

Detail™

Pre-Impression Cleansing Gel

10% EDTA Cleanser

A truly clean, uncontaminated prep site will produce a more accurate impression. Too often, hemostatic agents, coagulum or residual monomer from provisional fabrication can interfere with impression setting or result in less than ideal impression accuracy.

Detail™ Pre-Impression Cleansing Gel guarantees a clean, decontaminated prep site and consequently, an accurate impression. In just 7-10 seconds, **Detail** will lift organic debris and contaminants like sulfur and methacrylate residue from your prep site. Rinse away **Detail**, and the impression surface left behind will be squeaky-clean.

Detail does more than just clean the prep site, however. **Detail** also lowers the surface tension of the preparation so your impression material contacts the surface more intimately. The result? The most detailed, bubble-free and accurate impression possible.

For the best and most efficient impression procedure, use **Detail** in conjunction with 3rd Generation **Affinity™** Hydroactive Impression Material.



Fig. 1 – Contaminated preparation site - preventing an accurate and properly cured impression.





Fig. 2 – Clean preparation site. In just 7-10 seconds, Detail removes contaminants and surface debris, leaving the dentin tubule plugs in place, and allowing for a clean, decontaminated preparation site.


CLINICAL TECHNIQUE

Dentistry courtesy of Dr. Brad Carson. Photography courtesy of Dr. Gildo Santos.



In only 7-10 seconds, Detail leaves behind not only a clean preparation, but also a hydroactive surface – both will yield the cleanest, most detailed impression you've ever seen!

Dentistry and photography courtesy of Dr. Steve Tsotsos.



Prep site covered with Cerec powder (titanium dioxide opaquing medium) for Cerec imaging. It is critical to remove the powder for good adhesion.

Prep site after a 10 second application of Detail.

DETAIL 4-PACK

Contains: 4 x 3mL syringes Detail, 50 x flocked needle tips, instructions/SDS

TISSUE GOO/DETAIL TRIAL PACK

Contains: 2 x 3mL syringes Tissue Goo, 2 x 3mL syringes Detail, 25 x 20 gauge needle tips, instructions/SDS