


Introducing True™ Dental Dam



Isolate With Confidence.

True™ Dental Dam provides the fit, seal and isolation you demand for better clinical outcomes.



*Dentistry and photography courtesy of
Bob Margeas, DDS*

Simple Placement

True Dental Dam is easily punched, stretched and placed supporting effective and efficient isolation set-up procedures. Experience how smoothly **True Dental Dams** glide through contacts and tuck into the sulcus allowing for true placement confidence.

Durability & Protection

Exceptionally durable and tear resistant, **True Dental Dams** deliver the seal and tissue protection you expect and help to reduce contamination to create a clean, clear operative field.

Available in medium gauge, 5" x 5" or 6" x 6" latex (blue) and latex-free (purple) options to provide better color contrast with dentition.

Physical Properties

True Dental Dams - Latex-Free

Thickness	0.22mm
Tensile Strength	17.5 MPa
Elongation at Break	1100%
Tear Strength	55 N/mm

True Dental Dams - Latex

Thickness	0.23mm
Tensile Strength	23.5 MPa
Elongation at Break	700%
Tear Strength	85 N/mm



TRUE™ DENTAL DAM - LATEX
Medium gauge
5" x 5" 52-pack (**838252**)
6" x 6" 36-pack (**838236**)

NEW

TRUE™ DENTAL DAM - LATEX FREE
Medium gauge
5" x 5" 50-pack (**838200**)
6" x 6" 50-pack (**838201**)

NEW