


CounterFit™ II

Multi-Purpose Replication Silicone



CounterFit™ II is a silicone alginate replacement material, designed for all alginate-type impressions. Overcoming many clinical problems associated with alginate, such as: tearing, voids, lack of impression detail, and single pour restrictions; Counter Fit II:

- Flows beautifully to capture exquisite detail for all alginate type impressions.
- Exhibits an extremely high tear strength compared with traditional alginate and it cures to an ideal durometer to ensure interproximal tags do not fold during reinsertion.
- Once polymerized, can be poured more than once and reused, if necessary.
- Thixotropic to help prevent patient gagging.
- CounterFit II impressions can be stored as long as any VPS impression.




Applications

CounterFit II is ideally suited for any application requiring occlusal accuracy:

- Opposing arch impressions
- Bleaching tray & athletic mouthguard fabrication
- Study models
- Orthodontic models, appliances, and retainers
- Matrix material for long span provisional crown and bridge restorations
- TMJ splints
- All alginate-type impressions

MAXIMUM WORKING TIME:	1:10 min.
RECOMMENDED MINIMUM INTRAORAL SET TIME:	1:30 min
TOTAL SET TIME:	2:30 min.

Dentistry and photography courtesy of Len Boksmann, DDS, BSc, FADI, FICD



COUNTERFIT II 50ML 2-PACK (275602) Contains: 2 x 50mL Cartridges, 6 x 7.5 Mixing Tips, Instructions/SDS
COUNTERFIT II 50ML 12-PACK (275612) Contains: 12 x 50mL Cartridges, Instructions/SDS *mixing tips not included
COUNTERFIT II 50ML 36-PACK (275636) Contains: 36 x 50mL Cartridges, Instructions/SDS *mixing tips not included
COUNTERFIT II 380ML MEGAMIX CARTRIDGE (275382) Contains: 1 x 380mL Cartridge, Instructions/SDS *mixing tips not included
COUNTERFIT II 120ML 2-PACK (275603) Contains: 2 x 120mL Cartridges, 6 Dynamic Mixing Tips, Instructions/SDS

To view the latest news, education and more, visit:

www.clinicianschoice.com

