

Glisten™

Provisional Resin Glaze



Glisten™ provisional polish minimizes your polishing time while providing a beautiful, glistening shine on your provisional veneers, crowns, bridges, inlays and onlays.

- Once cured gives the provisionals a high gloss shine, like a natural tooth.
- Can be applied and cured extraorally or intraorally due to its resin-based formula.
- Provides instant polish for all provisional restorations.

CLINICAL TIP:

If GLISTEN is placed prior to cementation, it will prevent excess cement from adhering to the outside of the crown - saving finishing time.

CLINICAL TECHNIQUE


Intraoral Clinical Technique

Dentistry and photography courtesy of Len Boksmann, DDS, BSc, FADI, FICD and Brad Carson, DDS



Extraoral Clinical Technique

Dentistry and photography courtesy of Gildo Santos, DDS



GLISTEN (151207)
1 x 5mL bottle, instructions/SDS

To view the latest news, education and more, visit:

www.clinicianschoice.com   

