

Dr. Ronald Jordan

Composite Instrument Series



Comfortable and Reliable

The Dr. Ronald Jordan (REJ) Composite Instrument Series allows for enhanced handling, no material pull back, and perfect placement. Made with proprietary stainless steel and coated in non-porous Titanium Nitride, the instrument tips provide strength, flexibility, and a non-stick, smooth surface that enables accurate composite placement.

- Series of 8 non-stick composite instruments inspired by Dr. Ronald Jordan, a pioneer of esthetic dentistry.
- REJ Instruments feature large, autoclavable, comfortable Radel handles for precise placement and contouring of composite material.
- Titanium Nitride tips are non-stick for ease of placement.

Anterior Instruments:



REJ #01
Double paddle



REJ #02
Recommended fir initial labial compaction and sculpting. Blades are slightly longer and wider than REJ#01.



REJ #03
Thin, flexible blade reaches under margins to hone and sculpt. Ideal for working interproximally and packing retraction cord.



REJ #04
Designed for interproximal (not for packing composite). Ideal for shaping composite and embrasure areas of diastema closure restorations.

Anterior/Posterior Instrument:



REJ #10
Paddle with condenser ends with a blade slightly longer and wider than REJ#02.



REJ #20
Double-ended condenser.



REJ #21
Double-ended condenser, larger than REJ#20.



REJ #22
The rounded end is designed for use on occlusal anatomy and delicate shaping of the inclines.

Posterior Instruments:

Rounded tips spread composite material better and prevent pull back from sharp internal line angles.



PLEASE NOTE:
Once instruments are opened, they cannot be returned.

**DR. RONALD JORDAN
COMPOSITE INSTRUMENT SERIES**
REJ #01 (06101)
REJ #02 (06102)
REJ #03 (06103)
REJ #04 (06104)
REJ #10 (06110)
REJ #20 (06120)
REJ #21 (06121)
REJ #22 (06122)

To view the latest news, education and more, visit:

www.clinicianschoice.com   

 **Clinician's Choice®**
Designed for better dentistry.