

# Re-Cord™

## Knitted Retraction Cord








**Re-Cord™** non-impregnated retraction cord is unique, 100% cotton knitted cord that is designed to compress easily when packed into the sulcus.

- Following placement, the knitted loops of the retraction cord open which exerts a gentle continuous outward force, minimizing the trauma to the soft tissue, preventing gingival recession and biological width impingement.
- The unique design allows it to carry significantly greater quantities of hemostatic solution as well as provide an excellent medium for hemostatic gel absorption and placement.
- No memory facilitates easy packing and stays put when packed into the sulcus.
- Available in 6 color-coded sizes, each bottle contains 8ft/244cm of cord.

## CLINICAL TECHNIQUE

*Dentistry and photography courtesy of Brad Carson, DDS.*



### RE-CORD KNITTED RETRACTION CORD

#000 (649101)  
#00 (649102)  
#0 (649103)  
#1 (649104)  
#2 (649105)  
#3 (649106)



To view the latest news, education and more, visit:

[www.clinicianschoice.com](http://www.clinicianschoice.com)   

 **Clinician's Choice®**  
Designed for better dentistry.