

# PowerMix™

## Automatic Impression Material Dispenser

### IMPRESSIONS

#### DISPENSING POWER AT YOUR FINGERTIPS

Fast, effortless, and portable, the all-new **PowerMix™ Automatic VPS Dispenser** from **CLINICIAN'S CHOICE**, delivers perfectly mixed, void-free impression material with the press of a button. Now you can fill a full arch impression tray in 10 seconds or less with **PowerMix** - no more hand pumping.

**Fast:** Provides an even mix, eliminating voids, and in a fraction of the time (10 seconds).

**Convenient:** Effortless tray loading without hand fatigue.

**Portable:** Easily transportable from one operatory to the next.

**Lightweight:** Weighs 2lbs with full cartridge.

**Versatile:** Extrude material chairside or dispense via its small stationary base.

**Rechargeable:** A single battery supports 8 fully loaded 120mL cartridges before needing a re-charge.

#### COMBINE POWERMIX WITH 3<sup>RD</sup> GENERATION AFFINITY™ IMPRESSION MATERIAL

Clinically proven for ultra-detailed accuracy, **AFFINITY Heavy Body** is designed with an ideal flow and durometer for accurate full arch tray impressions, while **AFFINITY InFlex**, an ultra rigid tray material is designed specifically to prevent distortion with the dual-arch technique. Both **AFFINITY** tray materials are available in 120mL PowerMix cartridges, which are designed to reduce plastic waste and impression material left in the cartridge. Each 120mL **PowerMix** cartridge is enough for 3 full arch trays or 6-10 quadrant impression trays. Together, impression taking has never been more efficient or more accurate.



#### WORKING & SET TIMES

##### AVAILABLE WORKING TIME

Regular Set 1:45 min  
Fast Set 1:10 min

##### MINIMUM INTRAORAL SET TIME

Regular Set 2:30 min  
Fast Set 1:45 min

2:30 min reflects the minimum intraoral set time. Our recommendation is to allow 3:00 min for a full intraoral set.

1:45 min reflects the minimum intraoral set time. Our recommendation is 2:00 min for full intraoral set.

#### POWERMIX AUTOMATIC DISPENSER

PowerMix Dispenser Unit, Stand, AC Adapter, Instructions

#### AFFINITY INFLEX REGULAR SET 120mL 2-PACK

Contains: 2 x 120mL Cartridges, 6 Dynamic Mixing Tips, Instructions/SDS

#### AFFINITY INFLEX FAST SET 120mL 2-PACK

Contains: 2 x 120mL Cartridges, 6 Dynamic Mixing Tips, Instructions/SDS

#### AFFINITY HEAVY BODY REGULAR SET 120mL 2-PACK

Contains: 2 x 120mL Cartridges, 6 Dynamic Mixing Tips, Instructions/SDS

#### AFFINITY HEAVY BODY FAST SET 120mL 2-PACK

Contains: 2 x 120mL Cartridges, 6 Dynamic Mixing Tips, Instructions/SDS

#### DYNAMIC MIXING TIP FOR POWERMIX

50 x Dynamic Mixing Tips

**CLINICIAN'S CHOICE®**  
DENTAL PRODUCTS INC.