




## 5 REASONS TO BUY

**1** Diamond-impregnated, flexible polishing spirals are ideal for use on all composite, lithium disilicate, bisacryl, and most ceramic and zirconia materials.

**5** Autoclavable up to 30 times, A.S.A.P. All Surface Access Polishers can withstand wear while providing an efficient and consistent polish use after use.



**2** Designed for optimal adaptation to any surface, including difficult to reach areas.

**3** Simple, two-step system consistently provides a high luster and shine in less than 60 seconds.

**4** The A.S.A.P. Pre-polisher gently diminishes small surface defects without affecting anatomy, and then the A.S.A.P. Final High Shine Polisher provides a fast, lifelike polish.

## A.S.A.P.® All Surface Access Polishers



After trying the A.S.A.P. All Surface Access Polishers, they instantly became a staple product in our operatory, helping us achieve a remarkable shine on all of our restorations. The combination of using the A.S.A.P. Pre-polisher followed by the Final High Shine Polisher creates a unique polish on composites that rivals ceramics.

Sunny Viridi, DMD  
Winnipeg, Manitoba, Canada



[clinicianschoice.com/product/asap-all-surface-access-polishers/](http://clinicianschoice.com/product/asap-all-surface-access-polishers/)

### SPECIFICATIONS

Spiral wheel shape creates an amazing shine in less than 1 minute.

Polishers offered in two sizes: small (10 mm) and large (14 mm).

Autoclavable up to 30 times.



Clinician's Choice  
clinicianschoice.com • 800-265-3444