

# A.S.A.P.<sup>®</sup>

## All Surface Access Polishers



### More Shine. Less Time.

Diamond-impregnated **A.S.A.P.<sup>®</sup> Polishers** provide all surface access polishing and are ideal for use on all composite surfaces, lithium disilicate, bisacryls, most ceramics, and zirconia.

- The A.S.A.P. Pre-polisher (purple) is embedded with diamond particles (44 microns) and gently diminishes small surface defects without affecting anatomy, and prepares the surface for a final high gloss shine.
- The A.S.A.P. Final High Shine Polisher (peach) is embedded with small, 3-6 micron sized diamond particles to provide a life-like polish in as little as 20 seconds.
- Autoclavable up to 30 times, A.S.A.P. Polishers can withstand wear while providing an efficient and consistent polish use after use.


### All Surface Access Polishing



### Proven Performance

A recent university-based study,\*\* comparing the gloss produced by A.S.A.P. Polishers (Clinician's Choice<sup>®</sup>) and Enhance<sup>®</sup>/PoGo<sup>®</sup> (Dentsply Sirona) found that A.S.A.P. Polishers produced a higher gloss results, after 30 seconds, than Enhance and Pogo did after 90 seconds.

### Maximum Gloss Per System




Lawson, N. University of Alabama at Birmingham, School of Dentistry. Gloss of composites with the use of several polishers. April 2022.

\*Enhance<sup>®</sup> and PoGo<sup>®</sup> are registered trademarks of Dentsply Sirona. \*\* Data on file.

“A.S.A.P. Polishers save chair time by quickly polishing porcelain and composite restorations in **less than a minute.**”

- Danièle Larose, DMD, AAACD





**CLINICAL Instant Replay**

**In a rush?**  
Watch this quick technique tip!

POWERED BY  
**Hands-OnLive**

Scan to watch

**Achieving a Superior Polish with A.S.A.P. All Surface Access Polishers**

Bob Margeas, DDS

**CLINICAL Instant Replay**

**In a rush?**  
Watch this quick technique tip!

POWERED BY  
**Hands-OnLive**

Scan to watch

**Composite Finishing Technique for an Esthetic Result**

Bob Margeas, DDS



**Request a complimentary sample:**  
[sample.clinicianchoice.com/asap-polisher-family](http://sample.clinicianchoice.com/asap-polisher-family)



**Scan to Read**  
A.S.A.P. Polishers Dental Advisor Evaluation



<p><b>A.S.A.P. STARTER KIT (130917)</b> Contains: 3 x Pre-polishers Small, 3 x Pre-polishers Large, 3 x Final High Shine Polishers Small, 3 x Final High Shine Polishers Large</p>	<p><b>A.S.A.P. TRIAL 2-PACK</b> Contains: 1 x Pre-polisher, 1 x Final High Shine Polisher Small (130910) Large (130911)</p>
<p><b>A.S.A.P. REFILL 6-PACK</b> Contains: 3 x Pre-polishers, 3 x Final High Shine Polishers Small (130913) Large (130912)</p>	<p><b>A.S.A.P. BULK PACK SMALL PRE-POLISHER (130927)</b> Contains: 12 x Pre-polishers Small</p>
<p><b>A.S.A.P. PRE-POLISHER 3-PACK</b> Small (131011) Large (131012)</p>	<p><b>A.S.A.P. BULK PACK LARGE PRE-POLISHER (130928)</b> Contains: 12 x Pre-polishers Large</p>
<p><b>A.S.A.P. FINAL HIGH SHINE POLISHER 3-PACK</b> Small (132013) Large (132014)</p>	<p><b>A.S.A.P. BULK PACK SMALL FINAL HIGH SHINE POLISHER (130925)</b> Contains: 12 x Final High Shine Polishers Small</p>
	<p><b>A.S.A.P. BULK PACK LARGE FINAL HIGH SHINE POLISHER (130926)</b> Contains: 12 x Final High Shine Polishers Large</p>

To view the latest news, education and more, visit:

[www.clinicianchoice.com](http://www.clinicianchoice.com)