



# A.S.A.P.®

All Surface Access Polishers



## More Shine. Less Time.

A.S.A.P.® Polishers are a simple, two-step polishing system that provides a high luster and shine in as little as 30 seconds.

- **All Surface Access:** The spiral shape allows access to all surfaces of the restoration; simplifying the procedure and reducing the need to inventory multiple systems.
- **Fast, High Shine:** Both polishers in the A.S.A.P. system are diamond-impregnated to provide a 67% faster polishing time and an 86% higher gloss than a leading polishing system<sup>1</sup>.
- **Durable:** A.S.A.P. Polishers can be autoclaved up to 30 cycles, making them remarkably durable and cost-effective\*


\*Placing excessive pressure on the polisher will lead to premature degradation. The polishers will also wear faster on harder indirect surfaces.

<sup>1</sup> University of Alabama at Birmingham, School of Dentistry April 2022. Gloss of composites with the use of two polishing systems.

## Proven Performance

A recent university-based study,\*\* comparing the gloss produced by A.S.A.P. Polishers (Clinician's Choice®) and Enhance®/PoGo® (Dentsply Sirona) found that A.S.A.P. Polishers produced a higher gloss results, after 30 seconds, than Enhance and Pogo did after 90 seconds.

### Maximum Gloss Per System



“A.S.A.P. Polishers save chair time by quickly polishing porcelain and composite restorations in **less than a minute.**”

- Danièle Larose, DMD, AAACD





**CLINICAL Instant Replay**

**In a rush?**  
Watch this quick technique tip!

POWERED BY  
**Hands-OnLive**

Scan to watch

**Achieving a Superior Polish with A.S.A.P. All Surface Access Polishers**

Bob Margeas, DDS

**CLINICAL Instant Replay**

**In a rush?**  
Watch this quick technique tip!

POWERED BY  
**Hands-OnLive**

Scan to watch

**Composite Finishing Technique for an Esthetic Result**

Bob Margeas, DDS



**Request a complimentary sample:**  
[sample.clinicianchoice.com/asap-polisher-family](http://sample.clinicianchoice.com/asap-polisher-family)



**Scan to Read**  
Dental Advisor Editor's Choice ++++½ Evaluation



**A.S.A.P. STARTER KIT (130917)**  
Contains: 3 x Pre-polishers Small, 3 x Pre-polishers Large, 3 x Final High Shine Polishers Small, 3 x Final High Shine Polishers Large

**A.S.A.P. REFILL 6-PACK**  
Contains: 3 x Pre-polishers, 3 x Final High Shine Polishers  
Small (130913)  
Large (130912)

**A.S.A.P. PRE-POLISHER 3-PACK**  
Small (131011)  
Large (131012)

**A.S.A.P. FINAL HIGH SHINE POLISHER 3-PACK**  
Small (132013)  
Large (132014)

**A.S.A.P. TRIAL 2-PACK**  
Contains: 1 x Pre-polisher, 1 x Final High Shine Polisher  
Small (130910)  
Large (130911)

**A.S.A.P. BULK PACK SMALL PRE-POLISHER (130927)**  
Contains: 12 x Pre-polishers Small

**A.S.A.P. BULK PACK LARGE PRE-POLISHER (130928)**  
Contains: 12 x Pre-polishers Large

**A.S.A.P. BULK PACK SMALL FINAL HIGH SHINE POLISHER (130925)**  
Contains: 12 x Final High Shine Polishers Small

**A.S.A.P. BULK PACK LARGE FINAL HIGH SHINE POLISHER (130926)**  
Contains: 12 x Final High Shine Polishers Large

To view the latest news, education and more, visit:

[www.clinicianchoice.com](http://www.clinicianchoice.com)