



True™ Dental Dam



Isolate With Confidence.

True™ Dental Dam's ideal surface texture and elasticity allows it to glide smoothly through contacts and tuck into the sulcus for easy placement and optimal marginal seal.

- **Tear Resistant:** Outstanding durability and tear resistance for optimized treatment times.
- **Effortless Placement:** Offered in two pre-cut sizes in medium and heavy gauges to provide ease of placement for any procedure.
- **Isolation and Comfort:** With latex and latex free options, True Dental Dam improves infection control, visibility, and maximizes patient comfort.

Physical Properties	Medium Gauge		Heavy Gauge	
	LATEX	LATEX FREE	LATEX	LATEX FREE
THICKNESS	0.23mm	0.22mm	0.26mm	0.26mm
TENSILE STRENGTH	23.5 MPa	17.5 MPa	23.5 MPa	17.5 MPa
ELONGATION AT BREAK	700%	1100%	700%	1100%
TEAR STRENGTH	85 N/mm	55 N/mm	85 N/mm	55 N/mm

“The True Dental dam is **easy to place** and allows for an **excellent seal and isolation.**”

- Bob Margeas, DDS



Request a complimentary sample:
sample.clinicianchoice.com/true-dental-dam



Scan to Read
 Dental Advisor Editor's Choice and +++++ Evaluation



LATEX (BLUE) - MEDIUM GAUGE
 5" x 5" 52-pack (838252)
 6" x 6" 36-pack (838236)

LATEX (BLUE) - HEAVY GAUGE
 5" x 5" 52-pack (838253)
 6" x 6" 36-pack (838237)

LATEX FREE (PURPLE) - MEDIUM GAUGE
 5" x 5" 50-pack (838200)
 6" x 6" 50-pack (838201)

LATEX FREE (PURPLE) - HEAVY GAUGE
 5" x 5" 50-pack (838202)
 6" x 6" 50-pack (838203)

To view the latest news, education and more, visit:

www.clinicianschoice.com   

 **Clinician's Choice®**
 Designed for better dentistry.