

ROUND

				
	SHAPE			
	MAJOR DIAMETER (mm)	1.4	1.6	2.1
	HEAD LENGTH (mm)	1.4	1.6	2.1

AVAILABLE IN: COARSE MEDIUM

ROUND END CYLINDER

			
	SHAPE		
	MAJOR DIAMETER (mm)	1.2	1.4
	HEAD LENGTH (mm)	7.0	8.0


AVAILABLE IN: COARSE

ROUND END TAPER

												
	SHAPE											
	MAJOR DIAMETER (mm)	1.2	1.4	1.6	1.8	1.8	2.1	2.5	1.6	1.8	1.8	2.1
	HEAD LENGTH (mm)	8.0	8.0	8.0	8.0	8.0	8.0	8.0	9.0	9.0	10.0	9.0

AVAILABLE IN: COARSE MEDIUM FINE

BEVEL END CYLINDER

		
	SHAPE	
	MAJOR DIAMETER (mm)	1.2
	HEAD LENGTH (mm)	8.0

AVAILABLE IN: COARSE MEDIUM FINE

FOOTBALL

			
	SHAPE		
	MAJOR DIAMETER (mm)	1.8	2.3
	HEAD LENGTH (mm)	3.4	4.3

AVAILABLE IN: COARSE MEDIUM FINE

POINTED FOOTBALL

				
				
	SHAPE	368-018	368-023	US 285.5
	MAJOR DIAMETER (mm)	1.8	2.3	2.3
	HEAD LENGTH (mm)	4.6	5.0	5.0

AVAILABLE IN: COARSE MEDIUM FINE

BARREL

		
	SHAPE	
	MAJOR DIAMETER (mm)	3.3
	HEAD LENGTH (mm)	5.5

AVAILABLE IN: COARSE MEDIUM

MODIFIED TAPER

			
	SHAPE		
	MAJOR DIAMETER (mm)	1.6	1.6
	HEAD LENGTH (mm)	8.0	8.0

AVAILABLE IN: COARSE MEDIUM FINE

FLAT END TAPER

		
	SHAPE	
	MAJOR DIAMETER (mm)	1.8
	HEAD LENGTH (mm)	8.0

AVAILABLE IN: COARSE MEDIUM FINE

FLAME

			
	SHAPE		
	MAJOR DIAMETER (mm)	1.4	1.4
	HEAD LENGTH (mm)	8.0	8.0

AVAILABLE IN: COARSE MEDIUM FINE

NEEDLE

		
	SHAPE	
	MAJOR DIAMETER (mm)	1.4
	HEAD LENGTH (mm)	8.0

AVAILABLE IN: COARSE