




Evanesce Silk®

Nano-Enhanced Universal Restorative

Strong. Simple. Beautiful.
Smooth, Sculptable Handling.

Evanesce Silk® is a nano-enhanced universal restorative composite that features silky smooth handling, available in high demand shades and opacities that are ideal for the most common restorative cases.



Dentistry and photography courtesy of Fred Peck, DDS, FAACD

“Evanesce Silk handles well without sticking or slumping, and provides a stunning polish in under 90 seconds.”

- Fred Peck, DDS, FAACD



Silky Smooth Handling

Evanesce Silk, with a softer consistency than Evanesce®, is easy to place and manipulate yet will not slump, providing optimal procedure efficiency.

Fast, High Gloss Finish

The unique blend of Evanesce fillers results in a fast, high gloss polish that rivals porcelain in as little as one minute.



Dentistry and photography courtesy of David Chan, DMD, AAACD


Exceptional Esthetics

Evanesce Silk blends seamlessly with natural tooth structure and offers 4 enamel, 8 universal, and 2 dentin shade/opacity options.



	ENAMEL (80% Opacity)	UNIVERSAL (85% Opacity)	DENTIN (90% Opacity)
A1	●	●	
A2	●	●	●
A3	●	●	●
A3.5		●	
B1	●	●	
B2		●	
C2		●	
BL1		●	

Available in ● 0.3 g Single Dose 10-pack

● 0.3 g Single Dose 20-pack



Technical Profile



Evanesce Silk	
Filler Content by Weight	81%
Filler Content by Volume	66%
Shrinkage in Water (<24 hrs)	2.3%
Flexural Strength (MPa)	140 MPa
Elasticity Modulus (GPa)	10.4 GPa
Compressive Strength (MPa)	340 MPa
Radiopacity (% Al.)	200% Al.
Surface Roughness (μm)	0.07μm

Related Products



MPa Universal™

Light-Cure Adhesive

MPa Universal™ is a single-component, light-cure adhesive developed for use with self-etch, selective etch or total etch bonding techniques.



A.S.A.P.®

All Surface Access Polishers

A.S.A.P.® Polishers are a simple and effective, two-step intraoral polishing system that provides a high luster and shine in as little as 30 seconds.



Scan to learn more
about Evanesce Silk



EVANESCE SILK SINGLE DOSE 10PK

Contains: 10 x 0.3 g single dose,
instructions/SDS

A1E (323570)
A2E (323571)
A3E (323572)
B1E (323573)
A2D (323574)
A3D (323575)
B2U (323576)

EVANESCE SILK SINGLE DOSE 20PK

Contains: 20 x 0.3 g single dose,
instructions/SDS

A1U (323563)
A2U (323564)
A3U (323562)
A3.5U (323565)
B1U (323567)
C2U (323577)
BL1U (323569)

*Coming Soon

To view the latest news, education and more, visit:
www.clinicianschoice.com



 **Clinician's Choice®**
Designed for better dentistry.