

# Evanesce® Injectable

Injectable Restorative

## Versatile. Esthetic. Strong.

Introducing **Evanesce® Injectable** — the next-generation flowable composite ideal for the injection molding technique, providing excellent handling and esthetic properties.



Dentistry and photography courtesy of Jennifer Bell, DDS, MAGD, FICD.


“Evanesce Injectable delivers seamless blending, easy handling, and beautiful esthetics—it’s my go-to composite for natural anterior restorations.”

- Jennifer Bell, DDS, MAGD, FICD




## Versatile Applications

With exceptional stackability, control, and adaptation, Evanesce Injectable is ideal for the injection molding technique as well as traditional flowable indications.



## Reliable Strength


Achieves greater fracture toughness than the leading injectable composite<sup>1</sup>, while maintaining non-slumping, stackable handling properties required for injection molding placement.



<sup>1</sup>Lawson, N. University of Alabama at Birmingham, School of Dentistry. Fracture Toughness of Injectable Composites. April 2024.

## Excellent Esthetics

Offered in 6 shades, including 4 that are designed to match the VITA classical A1-D4® shade guide, Evanesce Injectable makes shade selection easy.




 **Clinician's Choice®**

## Fast, High Gloss Polish

All Evanesce Injectable shades provide a fast, high gloss polish and mimic the appearance of natural dentition.

Gloss Progression of Dental Composites



## Technical Profile

Evanesce Injectable	
Filler Content by Weight	59%
Filler Content by Volume	46%
Shrinkage in Water (<24 hrs)	4.0%
Flexural Strength (MPa)	130 MPa
Elasticity Modulus (GPa)	5 GPa
Compressive Strength (MPa)	320 MPa
Radiopacity (% Al.)	150% Al.

## Procedure Related Products



### Affinity™ Crystal

Clear Matrix Impression Material

Affinity™ Crystal is a clear silicone matrix material designed for direct composite veneers, inlays, onlays, full arch restorations and provisional restorations using the injection molding technique.



### A.S.A.P.®

All Surface Access Polishers

A.S.A.P.® Polishers are a simple and effective, two-step intraoral polishing system that provides a high luster and shine in as little as 30 seconds.



Scan to learn more  
about Evanesce Injectable

To view the latest news, education and more, visit:

[www.clinicianschoice.com](http://www.clinicianschoice.com)



G-aenial™ is a trademark owned by GC America Inc. 3M® is a registered trademark owned by 3M Company. Filtek® is a registered trademark owned by Solvntum Intellectual Properties Company. EvoFlow® is a registered trademark owned by Ivoclar Vivadent, Inc. VITA classical A1-D4® shade guide is a registered trademark of VITA.

© 2025 Clinician's Choice. All rights reserved. 821300900 R:0 05/25SN



#### EVANESCE INJECTABLE

Contains: 2 x 1.6 g syringes, 10 mixing tips, instructions/SDS

A1 Universal (323526)  
A2 Universal (323527)  
A3 Universal (323530)  
B1 Universal (323528)  
BL1 Universal (323529)  
Enamel White (323516)

 **Clinician's Choice®**