

# Evanesce® Flow

Flowable Restorative

## Exceptional Control. Beautiful Esthetics.

Offering preferred handling, high strength, and amazing esthetics, **Evanesce® Flow** is your solution to achieving better clinical outcomes for all minimally-invasive direct restorations.



Dentistry and photography courtesy of Marc Geissberger, DDS, MA, BS, CPT

“Evanesce Flow has excellent optical properties resulting in an ideal blend with natural tooth structure.”

- Marc Geissberger, DDS, MA, BS, CPT



## Preferred Handling

Thixotropic for excellent control, Evanesce Flow will not slump on a vertical surface making it suitable for Class V and temporary veneer restorations.



Dentistry and photography courtesy of Megan Shelton, DMD

## High Strength

Highly filled to give strength, Evanesce Flow can be used as a stand-alone restorative or as a base/liner.


## Amazing Esthetics

Evanesce Flow is available in 3 opacities: Enamel (80%), Universal (85%), and Dentin (90%) and 10 VITA shades for procedural versatility.

### Evanesce Flow Shades and Opacities

Enamel Shades (80% Opacity)					
A1	A2	A3		B1	
Universal Shades (85% Opacity)					
A1	A2	A3	A4	B1	BL1
Dentin Shades (90% Opacity)					
		A2			

## Technical Profile



Evanesce Flow	
Filler Content by Weight	62%
Average Filler Particle Size	0.5µm
Flexural Strength	125 MPa
Flexural Modulus	5600 MPa
Compressive Strength	280 MPa
Radiopacity	150% Al.


## Procedure Related Products



### MPa Universal™

Light-Cure Adhesive

MPa Universal™ is a single-component, light-cure adhesive developed for use with self-etch, selective etch or total etch bonding techniques.



### A.S.A.P.®

All Surface Access Polishers

A.S.A.P.® Polishers are a simple and effective, two-step intraoral polishing system that provides a high luster and shine in as little as 30 seconds.



Scan to learn more  
about Evanesce Flow



#### EVANESCE FLOW 1.7 G SYRINGE

Contains: 1 x 1.7 g syringe, 5 x 19 gauge tips, instructions/SDS

A1 Universal (323517)  
A2 Universal (323518)  
A3 Universal (323519)  
A4 Universal (323520)  
B1 Universal (323521)  
BL1 Universal (323525)  
A1 Enamel (323501)  
A2 Enamel (323502)  
A3 Enamel (323503)  
B1 Enamel (323505)  
A2 Dentin (323511)

To view the latest news, education and more, visit:

[www.clinicianschoice.com](http://www.clinicianschoice.com)    

© 2025 Clinician's Choice. All rights reserved. 821302100 R:0 07/25SN

 **Clinician's Choice®**